


# Electronic Balance:



## Standard



## Optional



Model	WT2003GH	WT3003GH	WT5003GH
Max Cap.	200g	300g	500g
Readability	1mg		
Pan Size(mm)	Φ80		
Display	LCD(White Back Light)		
I/O	RS232		
Power	220V/110V		
Unit	g, ct, oz, lb, x, %, pcs		
Outline Dimension(cm)	35×21×36		
Packing Size(cm)	42×32×44		
G.W.(kg)	6		
Main Function	Unit Transforms, Counting, Percentage, Printing		
Standard	RS232, Weight		
Optional	Pritner(dot-matrix printer, thermal printer), under hook		