



شركة أفلو
للأجهزة الطبية والمختبرية



YX930D Electric Suction Apparatus

Performance characteristics:

- 1 Large flow oil-free lubricating vacuum pump, high pumping rate, no oil mist pollution, pump body no need daily maintenance.
- 2 Large diameter liquid storage bottle, with a sealing ring stopper, easy to open and close the stopper, easy to remove the dirty liquid in the bottle.
- 3 Overflow protection device can prevent liquid from entering the pipe and pump.
- 4 Use transparent non-toxic PVC suction hose to observe the liquid in the tube.
- 5 Equipped with air filter can prevent negative pressure pump pollution.
- 6 The machine works smoothly, low noise and long service life.
- 7 Manual switch and foot switch are connected in parallel, optional.

YX930D Electric Suction Apparatus

The main technical parameters

- 1 Limit negative pressure value: $\geq 0.09\text{mpa}$ (680mmHg)
- 2 Negative pressure adjustment range: 0.02~0.09MPa(150~680mmHg)
- 3 Pumping rate: $\geq 40\text{L/min}$
- 4 Liquid storage bottle: 2500ml \times 2 (glass)
- 5 The noise: 60 db (A) or less
- 6 Power supply: AC220V 50Hz
- 7 Input power: 250VA
- 8 Outer packing size: 46.5cm \times 42.5cm \times 92cm
- 9 Gross/net weight: 22.5/20kg

0783 330 0730 البصرة 0771 608 6395 بغداد
0751 815 6446 الموصل 0751 815 6446 أربيل
0770 562 4669 العمارة 0770 681 5747 سليمانية
0790 238 8611 الأنبار 0770 135 0023 كركوك
0774 024 2223 واسط

Arasaat Al-Hindia Dist. 929 St.21, Villa 32
Baghdad-IRAQ

شركة أفلو للاجهزة الطبية و المختبرية

@aflo ltd

aflo ltd

www.aflo ltd.com

المميزات

- 1 كاتم للصوت
- 2 لا يحتاج الى صيانة
- 3 مصمم للعمل لفترات طويلة
- 4 مقيم من قبل وزارة الصحة العراقية ومورد لعشرات المتشفيات الحكومية